labmodul

STERILE BENCH CLASS II ● LAF BÄNK KLASSE II

MARS

MARS CLASS II BENCH

- Operator & Environment Safety.
- Ergonomic Design.
- Energy Efficiency.
- Noise Reduction.

Microbiological Safety Cabinet

The primary purpose of an MSC (Microbiological Safety Cabinet), BSC (Biological Safety Cabinet) or LAF bench, as the safety benches are called, is to act as a means of protecting the laboratory worker and the environment from pathogens. All air that is blown out is HEPA-filtered when it leaves the biosafety cabinet when harmful microorganisms have been removed.

Safety cabinets are used daily for hours, where in addition to the protection of users and test materials, the human design factors (ergonomics) of the work become increasingly important. This includes lowering the noise level (for a more comfortable working environment), a height-adjustable stand or chair (for optimized seating position), panoramic side windows (more light in the cabinet), angled front door (provides better seating position) and even light sources (better view of the cabinet) to improve working conditions.

LAF benches are divided into different categories depending on the work done in them, it is important that the safety cabinet meets the requirements set in your laboratory regarding the safety of the user, the product and the environment.

A Class I Safety Bench protects the user and the environment and is used, among other things, for powder medicine handling and organic solvents.

A Class II BSC Safety Cabinet protects users, the environment and the product, and is used, for example, for microbiological work, medicine and cell work.

A Class III LAF bench provides the same protection as a Class II but is designed to work with infectious microorganisms used in Biosafety Risk Group 4 (eg Ebola), offering maximum protection for the environment, product and user. Also called isolator.

Sterile bench is a safety bench with laminar flow (= LAF) where the quality of the air is HEPA-filtered (H14 filter) and the product is optimally protected. This type is used, for example, in the handling of sterile objects and in the electronics industry.

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Biological safety cabinets class 2 - Mars

The Mars is a series of dual filter Class 2 cabinets which incorporates the latest in laminator technology, energy- saving designs and HEPA-filtration, with options that give the best protection for you, your products and your laboratory.

The Mars safety cabinet sets the standards in quality, design, and innovation that comes from a heritage of over 50 years' experience. Inventors of the first electrically operated angled front screen in the early 1990's, our engineering team are constantly pushing forward the frontiers of energy saving design concepts and general development of our cabinets.

How does it work?

The air flow enters the work chamber of the biological safety cabinet via the front aperture and continues under the worktop and up the back plenum, where 70% is recirculated through the main HEPA filter, to provide the down flow protection and 30 % exits out through the HEPA filter to exhaust.

The vertical laminar flow recirculation provides operator protection by way of the inflow, product protection by means of the down flow and environmental protection by means of the filtered exhaust.

Mars Class 2 is the safety cabinet of choice when working with microbiological samples, tissue culture, viral work, non-pathogenic manipulations, IVF, research and sampling procedures.

Available in four different working width sizes 900 mm, 1200 mm, 1500 mm or 1800 mm
The Mars safety cabinets are being manufactured & tested in compliance with EN 12469.

All models are available with a wide range of high-quality customized options to suit your individual needs or requirements.

Leaving a "greener-footprint"

The Mars safety cabinets take full advantage of the latest innovative energy-efficient fans and optimized air flows, contributing to lower energy consumption and reduction of the noise level.

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
Ultra clean environment – Safety first

- 110 mm HEPA filters with efficiency at 99,999 % against 0,3 µm particles, ensuring a clean sterile work chamber environment.
- Unique laminator/diffuser technology giving turbulent-free air flow, protecting your samples against particle contamination.
- Down flow & exhaust fans are constantly monitored via the flow sensors to ensure safe operating conditions are maintained.
- Alarms, both acoustic and visual for unsafe airflow conditions or any interruption.
- Angled pre-filter for easy inspection and filter exchange.
- Sectional 300 mm work tops allowing loads up to 50 kg. per 300 mm, easy removal for cleaning and decontamination of all surfaces.
- Easy to maintain as all service functions are performed from the front of the cabinet, including changing of the HEPA filters, pcb's and sensors, and adjustments/monitoring of alarms and fan speed.

Ultra Clean Environment – Safety First

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Energy saving benefits

New Features of the range:
Now available with the new innovative EC fan motors and optimized air flow.

For a Mars 1200 this equates to an energy consumption of 92* Watts at factory settings and 30 Watts in ECO mode. This incremental innovation makes the Mars cabinet range some of the lowest energy consuming cabinets on the market!

The innovation of using several low energy EC fans in the cabinet allows for a less restricted construction and a reduction of annual operating costs, whilst at the same time allowing for the use of the 110 mm HEPA filter, equating to a 50 % longer filter life.

- Low energy consumption result in less heat transmission to the work chamber and to the laboratory contributing significantly to a reduction of overall energy costs.

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

Annual energy savings when comparing Mars 1200 to old equivalent cabinets in the industry

	Power consumption/W	Energy consumption/kWh**	Costs in EUR/year***	Costs/15 years****
Old generation	420	739	229	3.437
Mars 1200	92	162	50	753
Reduction of overall energy costs				78 % saving, EUR 2.684

- Different types of electrical outlets can be ordered.
- Sectional work tops can be provided in different sizes from 300mm to 1800mm and also in AISI 316 to enable the MARS to be configured to your exact practical requirements.
- Sink with a water tap mounted into a solid one-piece worktop.
- A marble stone can be inserted into the work top section for a balance to enable weighing applications to 6 decimals.
- Built-in LAF-LCD screen mounted on the rear wall or alternatively magnetically mounted for easy removal or positioning to suit the operator's convenience and comfort.

Numerous other options are available, ranging from different valve types, the interior in AISI 316 instead of AISI 304, heated-plate zone system, PIR sensor, Bunsen burner, LED light, mounting of microscope etc.

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


Connection of your Mars cabinet to external ducting system

- HEPA filters do not trap gases and toxic chemicals that may be used as an adjunct to microbiological studies. In these applications it is important to select an exhaust ducting route to lead the exhausted air safely out of the laboratory to an external ventilation system.
- It is possible to connect the cabinet to an external ducting system via an anti-blow back valve. The extract connection can be either referred to as a hard duct or a thimble duct, where the latter is also extracting air from the laboratory.
- The compact construction of the Mars Class 2 cabinet makes it possible for it to be used – ducted either way – with a working height of up to 1000 mm in a laboratory with a ceiling height of 2,5 meters.

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LABMODUL A/S